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- Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
- Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
- Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.
- Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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- International Federation of Clinical Chemistry, Scientific committee. J Clin Chem Clin Biochem 1980 18: 521-534.
- Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
- Bergmeyer HU, Herder M, Rej R. Approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 2. IFCC Method for aspartate aminotransferase. J Clin Chem Clin Biochem 1986;24:49.

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Inter Assay – Between Run			
Sample	Mean (U/L)	SD (U/L)	CV%
Control Serum 1	42	1.2	2.85
Control Serum 2	161	1.9	1.18
Control Serum 3	202	2.5	1.23

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

$$Y=0.98(X)-2.2 \text{ U/L}; r=0.99$$

	For In Vitro Diagnostics Use Only
	Lot Number
	Catalogue Number
	Storage Temperature
	Expiry Date (Year / Month)
	Warning, Read Enclosed Documents
	Instructions For Use
	Manufactured By