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Kit Contents:	BXC0222A
R1 BCG Reagent	4 x 125 ml

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	Reagent Blank	Standard / Sample
Standard / Sample	---	10 μl
Reagent	1000 μl	1000 μl

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Inter Assay – Between Run			
Sample	Mean (g/dl)	SD (g/dl)	CV%
Control Serum 1	1.8	0.02	1.11
Control Serum 2	3.7	0.04	1.08
Control Serum 3	6.3	0.18	2.85

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- Grant GH., et al Amino Acids and Protein. In Fundamentals of Clinical Chemistry, Tietz N.W. Editor, Third Edition, WB Saunders Company Philadelphia USA, 328-329 1987
- Young DS. Effects of drugs on clinical laboratory tests, AACCC Press 2000
- Doumas BT, Watson WA, Biggs HG. Albumin standards and the measurement of serum albumin with bromocresol green. Clin Chem Acta 1971
- Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
- Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
- Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.
- Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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