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Intra Assay – Within run			
Sample	Mean(mg/l)	SD(mg/l)	CV %
Control Serum 1	0.2	0.01	5.0
Control Serum 2	0.6	0.04	6.66
Control Serum 3	1.1	0.06	5.45

Inter Assay – Between Run			
Sample	Mean(mg/l)	SD(mg/l)	CV %
Control Serum 1	0.25	0.02	8.0
Control Serum 2	0.71	0.05	7.04
Control Serum 3	1.2	0.07	5.83

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1. Teuscher, A & Richterich, P. Schweiz Med Wschr. 1971;101:345 & 390
2. Tietz NW, Clinical guide to Laboratory Tests, 3<sup>rd</sup> Edition.
3. Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
4. Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
5. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
6. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
7. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.
8. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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