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EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.

- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.
- Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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- Assmann G. At. what levels of total low-or high-density lipoprotein cholesteol should diet/drug therapy be initiated European guidelines. Amer J Cardiol 1990;65:11F
- Assmann G. Schriewer H. Schmitz G et al Qualification of high density lipoprotein cholesterol by precipitation with phosphotungstic acid/Mg/Cl₂ Clin Chem 1983;29:2026-2030
- AVP Fettstoffwechselstörungen, Therapieempfehlungen 1, 1st 1996:2-16
- Bablok W. et al. General Regression Procedure for Method Transformation. J Clin Chem Biochem 1988; 26:783-790
- Burstein M. Scholnick HR. Morfix R. Rapid method for the isolation of lipoproteins from human serum by precipitation with polyanions. J Lipid Res 1970;11:583-595
- Cohn J.S., McNamara J.R., Schaefer E.J.. Lipoprotein Cholesterol Concentrations in the plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459
- Glick M.R., Ryder K.W., Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
- Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
- Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
- Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document

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