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Kit Contents:	BXC0120A
R1: Buffer Reagent	1x30 ml
R2: MDH Reagent	1x2 ml
R3: Enzymatic Reagent	1 Vial
R4: Standard controls	1x5ml 2X3 ml

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- Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
- Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
- Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.
- Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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- Marty et al. Clin.chem. 30/7 .1231 (1984)
- Milan A, Cont A, Gacia-Raso A, Garases F (1987). Determination of citrate in urine by ClinChem 33.1259-1260. Simple direct photometry
- Top s, Yucel D Determination of citrate in urine by simple direct photometry. Clin chem 34.1658 (1988)

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Intra Assay – Within Run			
Sample	Mean (mg/l)	SD (mg/l)	CV%
Control 1	450	4.2	0.93
Control 2	610	3.5	0.57
Control 3	960	3.5	0.36

Inter Assay – Between Run			
Sample	Mean (mg/l)	SD (mg/l)	CV%
Control 1	465	4.7	1.01
Control 2	602	4.1	0.68
Control 3	966	4.4	0.45

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
	For in Vitro Diagnostics Use Only
	Lot Number
	Catalogue Number
	Storage Temperature
	Expiry Date (Year / Month)
	Warning, Read Enclosed Documents
	Instructions for Use
	Manufactured By