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- Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
- Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
- Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.
- Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

|   |                                   |
|---|-----------------------------------|
|  | For In Vitro Diagnostics Use Only |
|  | Lot Number                        |
|  | Catalogue Number                  |
|  | Storage Temperature               |
|  | Expiry Date (Year / Month)        |
|  | Warning, Read Enclosed Documents  |
|  | Instructions For Use              |
|  | Manufactured By                   |

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- Assmann G. At. what levels of total low-or high-density lipoprotein cholesteol should diet/drug therapy be initiated European guidelines. Amer J Cardiol 1990;65:11F
- Assmann G. Schriewer H. Schmitz G et al Qualification of high density lipoprotein cholesterol by precipitation with phosphotungstic acid/Mg/Cl<sub>2</sub> Clin Chem 1983;29:2026-2030
- AVP Fettstoffwechselstörungen, Therapieempfehlungen 1, 1<sup>st</sup> 1996:2-16
- Bablok W. et al. General Regression Procedure for Method Transformation. J Clin Chem Biochem 1988; 26:783-790
- Burstein M. Scholnick HR. Morfix R. Rapid method for the isolation of lipoproteins from human serum by precipitation with polyanions. J Lipid Res 1970;11:583-595
- Cohn J.S., McNamara J.R., Schaefer E.J.. Lipoprotein Cholesterol Concentrations in the plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459
- Glick M.R., Ryder K.W., Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474

| Sample          | Inter Assay – Between Run |           |      |
|-----------------|---------------------------|-----------|------|
|                 | Mean(mg/dl)               | SD(mg/dl) | CV % |
| Control Serum 1 | 22                        | 0.57      | 2.5  |
| Control Serum 2 | 48                        | 0.45      | 0.93 |
| Control Serum 3 | 75                        | 0.9       | 1.2  |

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