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$$\text{Phosphorus (mg/dl)} \times 0.32 = \text{Phosphorus (mmol/l)}$$

$$\text{Phosphorus (mg/dl)} \times 3.06 = \text{Phosphate (mg/dl)}$$

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$$\text{urine phosphorus (mg/24h)} = \frac{\text{urine ph. (mg/dl)} \times \text{urine volume (ml)}}{100}$$

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Intra Assay – Within Run			
Sample	Mean (mg/dl)	SD (mg/dl)	CV%
Control Serum 1	2.0	0.02	1.0
Control Serum 2	4.2	0.05	1.19
Control Serum 3	6.0	0.07	1.16

Inter Assay – Between Run			
Sample	Mean (mg/dl)	SD (mg/dl)	CV%
Control Serum 1	2.2	0.05	2.27
Control Serum 2	4.35	0.08	1.83
Control Serum 3	6.2	0.09	1.45

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6. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.

7. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.

8. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

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1. Tietz NW. *Clinical Guide to Laboratory Tests* 3rd Edition
2. Henry RJ. *Clinical Chemistry Principles and Techniques*
3. Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
4. Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
5. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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	For In Vitro Diagnostics Use Only
	Lot Number
	Catalogue Number
	Storage Temperature
	Expiry Date (Year / Month)
	Warning, Read Enclosed Documents
	Instructions For Use
	Manufactured By