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9. Passing H., Bablok W., A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
10. Pisani T., GebSKI C.P. Leary E.T., et al. Accurate Direct
11. Determination of Low-Density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
12. Rifai N., Warnick G.R., Mc Namara J.R., Belcher J.D., Grinstead G.F., Frantz Jr I.D. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 1992;38:150-160.
13. 4. Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
14. 5. Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
15. 6. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.

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1. Armstrong V., Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Ärztl. Lab. 1985;31:325-330.
2. Bablok W. et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
3. Bachorik P.S., Ross J.W. National cholesterol education program recommendations for measurement of low-density lipoprotein cholesterol: executive summary. Clin Chem 1995;41:1414-1420.
4. Cohn J.S., Mc Namara J.R., Schaefer E.J.. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
5. Cremer F., Nagel D., Mann H. et al. Ten-year follow-up results from the Goettingen Risk, Incidence and Prevalence Study (GR(PS). 1. Risk factors for myocardial infarction in a cohort of 5790 men. Atherosclerosis 1997;129:221-230.
6. Cremer F. Seidel D. Lipoproteinanalytik: Methodische Empfehlungen. Friedewald W.F., Levy R.I., Frederickson D.S. Estimation of LDL-Cholesterol Concentration without Use of the Preparative Ultra-centrifuge. Clin Chem 1972;18:499-502.
7. DG Klinische Chemie. Mitteilungen 1990;21:215-232.
8. Naito H.K., Strong J.P., Scott M.G., Roheim P.S., Asztalos B.F., Zilversmit D.B., Srinivasan S.R., Berenson G.S., Wilson P.W.F., Scanu A.M., Malikow M.R., Atherogenesis: current topics on etiology and risk factors. Clin Chem 1995;41:132-133 No. 1.

Sample	Inter Assay – Between Run		
	Mean(mg/dl)	SD(mg/dl)	CV %
Control Serum 1	115	0.91	0.79
Control Serum 2	136	1.21	0.88
Control Serum 3	180	1.51	0.83

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Recommended (desirable)
< 130 mg/dl
Moderate Risk
130-159 mg/dl
Risk
≥ 160 mg/dl

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