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$$\text{Total Protein (mg/24h)} = \frac{\text{Total Protein (mg/dl)} \times \text{Urine Volume (ml)}}{100}$$

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Intra Assay – Within Run			
Sample	Mean (mg/dl)	SD (mg/dl)	CV%
Control 1	21.4	0.31	1.44
Control 2	80	1.22	1.52
Control 3	121.5	2.65	2.18

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8. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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Sample	Inter Assay – Between Run		
	Mean (mg/dl)	SD (mg/dl)	CV%
Control 1	22.6	0.75	3.31
Control 2	81.5	1.86	2.28
Control 3	123.6	3.21	2.59

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Biorexfars Urinary Protein/CSF Standard Cat. No BXS0172

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1. Tietz NW., Clinical Guide to Laboratory Tests 3rd Edition
2. Henry RJ Clinical Chemistry Principles & Techniques 2nd Edition 1974
3. Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
4. Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
5. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
6. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
7. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.

	For In Vitro Diagnostics Use Only
	Lot Number
	Catalogue Number
	Storage Temperature
	Expiry Date (Year / Month)
	Warning, Read Enclosed Documents
	Instructions For Use
	Manufactured By