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Inter Assay – Between Run			
Sample	Mean(ug/dl)	SD(ug/dl)	CV %
Control 1	73	1.34	1.83
Control 2	155	2.85	1.83
Control 3	260	3.4	1.30

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15. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.

16. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.

17. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

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1. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988; 26:783-790.
2. Berthelot MPE. Reper Chem Appl. 1859:282.
3. Da Fonseca-Wollheim F. Z Klin Chem Klin Biochem. 1973;1 1 :421, 426.
4. Dokumentation Roche Diagnostics.
5. Hohorst HJ. Biochem Z. 1956;328:507.
6. Kaplan LA, Pesce AJ. Clinical Chemistry Theory, Analysis and Correlation. St. Louis, Mo: CV Mosby Co; 1984:1231.
7. Kirsten E, et al. Biochem Z. 1963;337:312.
8. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983; 21:709-720.
9. Prellwitz W, et al. Med Welt. 1976;27:1277.
10. Siegel JM, Montgomery GA. Arch Biochem/Biophys. 1956;82:288.
11. Tietz NW. Clinical Guide to Laboratory Tests. 3<sup>rd</sup> ed. Philadelphia, PA: W.B. Saunders Co; 1995:44.
12. Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
13. Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
14. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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