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**biorex fars**  
**AMYLASE CNPG3**  
**PRODUCT CODE: BXC0562**

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Institute; 2003.

9. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. CLSI document EP17-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2012.
10. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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| Sample          | Inter Assay – Between Run |          |      |
|-----------------|---------------------------|----------|------|
|                 | Mean (U/L)                | SD (U/L) | CV%  |
| Control Serum 1 | 176                       | 2.5      | 1.42 |
| Control Serum 2 | 381                       | 3.2      | 0.83 |
| Control Serum 3 | 689                       | 5.6      | 0.81 |

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

**Biorex fars Calibration Serum Cat. No BXC0321A**

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- Young DS, Effects of drugs on clinical Laboratory tests. 4th ed AACC Press. Washington DC ; 3-43 to 3-47 1995.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
- Greiling H, Gressner AM ed. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3. Auflage. Stuttgart/New York: Schattauer Verlag, 1995.
- Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
- Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
- Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards

|   |                                   |
|---|-----------------------------------|
|    | For In Vitro Diagnostics Use Only |
|    | Lot Number                        |
|    | Catalogue Number                  |
|    | Storage Temperature               |
|    | Expiry Date (Year / Month)        |
|   | Warning, Read Enclosed Documents  |
|  | Instructions For Use              |
|  | Manufactured By                   |