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Intra Assay – Within Run			
Sample	Mean(mg/L)	SD(mg/L)	CV %
Control 1	24.4	0.38	1.55
Control 2	41.0	0.65	1.58
Control 3	88.0	0.86	0.97

Inter Assay – Between Run			
Sample	Mean(mg/L)	SD(mg/L)	CV %
Control 1	25.5	0.77	3.01
Control 2	42.3	1.02	2.41
Control 3	89.2	1.1	1.23

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8. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP07-A2. Wayne,PA: Clinical and Laboratory Standards Institute; 2005.

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1. Elving, LD et al (1989) clin Chem: 35:2; 308.
2. Bakker AJ (1988) 34:1:82
3. 4. Evaluation of Precision of Quantitative Measurement; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2014.
4. Evaluation of Stability of In Vitro Diagnostics Reagents; Approved Guideline. CLSI document EP25-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2009.
5. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne,PA: Clinical and Laboratory Standards Institute; 2013.
6. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. CLSI document EP06-A. Wayne,PA: Clinical and Laboratory Standards Institute; 2003.
7. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.

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

Biorexfars Microalbumin Calibrator series cat. No BXC0320A (5x1ml)

	For In Vitro Diagnostics Use Only
	Lot Number
	Catalogue Number
	Storage Temperature
	Expiry Date (Year / Month)
	Warning, Read Enclosed Documents
	Instructions For Use
	Manufactured By